

# Cordless Suspense Point

(Battery Operated)

Indigo Care has designed a Cordless Suspense Point that will suspend all the cordless devices in the room. This device deactivates all of the Indigo Care's cordless sensor mats and neck/wrist pendants temporarily so that the care-givers can monitor the residents and/or perform general housekeeping with ease, without setting off the nurse call system.

The Suspense/Cancel point is designed with 2 buttons: a white SUSPENSE and a white CANCEL. Needing only a Lithium Coin Battery (CR 2032) to power, it can be installed anywhere, but ideally next to the door entrance. Simply use the included screw kit or double sided tape if preferred (not included).

Programming the two buttons to the RRU (Indigo Care's Cordless Receiver) is easy but must be done separately for SUSPENSE and CANCEL. See our Cordless Programming Instructions for complete details.

When the care-givers need to enter the room and suspend the RTU (Indigo Care's Cordless Transmitters), simply press the SUSPENSE button. The button will light-up momentarily. This will suspend all RTU devices in the room for 5 minutes (factory default setting).

You will notice that once suspended, the LED light on the RRU will change from a green continuous light to a blinking amber light. This is an indication that the RTU devices are in suspension mode.

Upon approaching the end of the 5 minute suspension, the RRU will give you a 10 second buzz warning. If you require more time, simply press the SUSPENSE button again, and the 5 minute cycle starts over.

If you leave the room before the 5 minutes lapses, press the CANCEL button upon your departure. Alternatively, it will automatically reset after 5 minutes.

At the end of the suspended time, the RRU will have a continuous green light, indicating that every RTU device in the room has been reset and is again active.

**Note:** The Suspense Point does *NOT* de-activate the corded pendant or other types of corded devices connected directly to the Nurse Call Point.



## Key Specifications:

- Australian Design
- C-Tick Approved
- Anti-Bacterial Silicon Buttons
- Raised, Easy to Press Buttons
- 2 Buttons: White SUSPENSE, White CANCEL
- ID Encoded RF Transmitter
- Momentary Back Light
- Input Supply: Lithium Coin Battery CR 2032
- Screw Kit included
- Optional Double Sided Tape
- Suits Standard Australian Electrical Boxes and Mounting Enclosures

Part number: MED-RTU-SUSP



HORIZON ELECTRONIC COMPONENTS Pty Ltd  
2/71 Rushdale Street, Knoxfield, VIC, 3180  
PH +61 3 9763-2988 FAX +61 3 9763-3811  
www.horizone.com.au  
ABN 36007302801

05/20



## ORDER INFORMATION

Fax or email your purchase orders to  
**03-9763-3811**  
Or email [sales@horizone.com.au](mailto:sales@horizone.com.au)