MULTI PURPOSE ALARM MANAGEMENT:

NORMALLY OPEN SWITCH

NORMALLY CLOSE SWITCH

MAGNETIC SWITCH

REDUCE THE RISK OF FALLS & INJURIES

PORTABLE

EASY TO USE

EASY TO MOVE FROM ROOM TO ROOM

INTEGRATES TO MOST NURSE CALL SYSTEMS

PLUG & PLAY OPERATION

DURABLE & RELIABLE

COST EFFECTIVE



Dry Sensor Monitor – Anti-Tamper



PRODUCT CODE: MED-NC-TS1011A (TGA Registration: 223902)

Indigo Care provides the Anti-Tamper Care Alarm Monitor that is designed to be used in conjunction with our floor mats and bed mats/pads. This is a multi-function fall management alarm which can be used as a monitor with magnetic switch, floor sensor mat and bed/chair sensor mat. The alarm will sound when a resident leaves their bed/chair or steps on the floor mat or when the magnetic switch is detached (when magnetic switch is in use). This care alarm can only be reset either by a Reset Key or an Infrared (IR) Remote Control. This eliminates the risk of the alarm being turned off by the resident. It can connect and alarm through your current nurse call system. There are two alarm sounds & a Low or High Volume setting. Alternatively, this alarm sound can be put on silent mode and used with your nurse call system.

KEY SPECIFICATIONS:

- Versatile mounting base & hanger included
- Tamper free reset & suspend by infrared remote control
- Suspend mode allows for patient repositioning with adjustable stand-by time
- Reset mode puts the alarm into stand-by until re-activated by next use
- 2 alarm sounds High & Low
- Silent mode for use with remote alarm signaling Nursecall systems
- 4 digits ID code to match remote control to alarm & eliminate interference
- Rubberized housing for protection against damage caused by dropping
- Versatile mounting base and hanger included
- Easy Clean Pull Cord
- Dual Power: Battery Operated 4 x AA (not included) or optional 6VDC Adaptor
- Warranty: 12 months

Horizon Electronic Components T/as INDIGO CARE 2/71 Rushdale Street, Knoxfield, VIC, 3180 PH +61 3 9763-2988 FAX +61 3 9763-3811 www.indigocare.com.au ABN 36007302801

ORDER INFORMATION
Fax or email your purchase orders to
03-9763-3811
Or email sales@indigocare.com.au