

Wireless Chair Mat

INTEGRATES TO MOST NURSE CALL SYSTEMS VIA OUR RECEIVERS

REDUCE THE RISK OF FALLS & INJURIES

SOLVENT WIPE-DISINFECTANT

LOW PROFILE, NON-SLIP MATERIAL

NO WIRES AROUND RESIDENT/PATIENT- NO TRIP HAZARDS

PORTABLE

EASY TO USE

EASY TO MOVE FROM ROOM TO ROOM

PLUG & PLAY OPERATION

DURABLE & RELIABLE

COST EFFECTIVE



TGA Registered

RCM Approved

AS/NZS 61000.6.3.3012

PRODUCT CODES: Chair Mat & Transmitter: MED-NC-C365X-WT-V1.1
Full Kit: Receiver Unit, Chair Mat & Transmitter: MED-NC-C365K-WT-V1.1

Indigo Care provides a pressure sensitive Wireless Chair Sensor Pad that is designed to immediately alert staff when a resident leaves their chair.

The transmitter is an exclusive ID encoded RF Transmitter which communicates with the Wireless Receiver Box (RRU). It allows the Chair sensor mat to wirelessly transmit an alarm to the call point. Simply connect the Transmitter to the chair sensor mat. Then teach in the Transmitter to RRU by following the instructions printed. The sensor mat will now trigger the call point when activated within the transmitting area. It is easily integrated to most nurse call systems. (Note: Transmitter and Bed mat can be ordered separately)

KEY SPECIFICATIONS:

Chair mat

- Dimensions: 380mm X 255mm
- Connection: RJ12 on 500mm curly cord
- Silver Printing, Resistance <1k Ohm
- White PVC, Thickness 0.20mm
- Warranty 12 months

Transmitter

- Connection: RJ 4P4C socket
- No External Antenna
- LED low battery indicator
- Transmits up to 30m (line of sight)
- Frequency: 433.92MHz
- Uses a 3V battery-CR2032
- Warranty: 12 months

PRODUCT ACCESSORIES:

- Wireless Receiver Box
- Wireless Floor Sensor Mat
- Wireless Bed Mat
- Wireless Waterproof pendant
- Wireless PIR Motion Sensor

Horizon Electronic Components

T/as INDIGO CARE

2/71 Rushdale Street, Knoxfield, VIC, 3180

PH +61 3 9763-2988 FAX +61 3 9763-3811

www.indigocare.com.au

ABN 36007302801



ORDER INFORMATION

Fax or email your purchase orders to
03-9763-3811

Or email sales@indigocare.com.au

