

SUITABLE FOR UNDER
THE MATTRESS

REDUCE THE RISK OF
FALLS & INJURIES

LOW PROFILE,
NON-SLIP MATERIAL

PORTABLE

EASY TO USE

EASY TO MOVE
FROM ROOM TO
ROOM

SOLVENT WIPE-
DISINFECTANT

INTEGRATES TO
MOST NURSE CALL
SYSTEMS

PLUG & PLAY
OPERATION

DURABLE & RELIABLE

COST EFFECTIVE

LONG PRESSURE PAD
(760mm X 110mm)



Under Mattress Sensor Pad



PRODUCT CODE: MED-B365-UM-4
(TGA Registration: 223902)

Indigo Care provides a pressure sensitive Under Mattress Sensor Pad that is designed for UNDER MATTRESS to immediately alert staff should a resident leaves the bed. This is most useful to residents who are at risk of falling out of their bed, or those who require assistance in arising from the bed. The sensor pad may not be suitable for some mattresses ex. air mattress. The bed sensor pad is connected to your existing nurse call system via the monitor and a Y-Cable. Once the resident removes themselves from the pad and tries to leave the bed, the sensor pad will activate the monitor's alarm and alert the staff. There are a set of instructions printed on the actual mat. Eyelets at the end of the pad are for easy storage.

Manufacturer/ Supplier recommendation:

Check the sensor pad is compatible for the mattress before purchased. Get a trial sample from Indigo Care.

KEY SPECIFICATIONS:

- **Dimensions: 760mm X 110mm**
- **Connection: RJ11 on 500mm curly cord**
- **Electrical Resistance <100 Ohm**
- **White PVC, Thickness 0.20mm**
- **Warranty: 12 months**

PRODUCT ACCESSORIES:

- Y Cable
- Anti-temper Monitor (use 4 x AA batteries-not included)
- AC Adaptor (optional)

Horizon Electronic Components T/as INDIGO CARE
2/71 Rushdale Street, Knoxfield, VIC, 3180
PH +61 3 9763-2988 FAX +61 3 9763-3811
www.indigocare.com.au
ABN 36007302801



ORDER INFORMATION

Fax or email your purchase orders to
03-9763-3811
Or email sales@indigocare.com.au